


Applications


- ◎ Analytical instruments
- ◎ Titration
- ◎ Filling
- ◎ Hydroponics
- ◎ Beauty devices



Pump Size



334K/BT



334K/TB

Features

- Convenient and quick tube installation, single pump head can hold 1 or 2 tubes to form a single or dual channel
- Can be connected in series with dual pump heads to double the flow rate
- Comes standard with 4-roller low-pulse rotor
- Tube compatible with 1.6mm wall thickness

Technical Specifications

Pump Model	334K/BT	334K/TB	
Motor Type	BT	Synchronous	
Power Supply	DC 24V	AC 220V,50Hz	
Power	≤30W	≤7W	
Net Weight	575g(Single Pump)/ 660g(Double Pump)	330g(Single Pump)/ 410g(Double Pump)	
Work Environment	Temperature range: 0℃ (no freezing) to 60℃, relative humidity <85%RH		
Output Pressure	≤0.17MPa		
Single pump, single channel - Flow Rate (unit:ml/min)			
rpm Tube#	300r	10r	30r
13#	10	0.3	1
14#	35	1.3	4
19#	80	2.5	8
16#	115	3.5	11.5
25#	210	7	22

