




FDA Certification
No.CANEC1925591702

GB 4806 Certification
No.CANEC2106280104

Features

- Various connector types available: straight head/extended straight head/ T-junction/hard tube connector/threaded connector/reducer straight head/reducer threaded connector/extended straight head/right angle connector
- Made of PP material, resistant to chemical contact
- Optional PVDF material for high corrosion resistance
- Temperature range: -13°C to 66°C, suitable for high-temperature sterilization
- Compatible with standard tubing, no limitations on external tubing material or hardness

Application



Optional connector specifications and compatible tubes

Y-shaped connector	Compatible tube	Reducer straight connector	Compatible tube	Threaded connector	Compatible tube
Φ1	1*1	Φ1.6-Φ2.4	2*1/14#-3*1/19#	Φ1.6-16#	2*1/3*1/16#/25#-ID4*1
Φ3.2	16#	Φ4.8-Φ3.2	25#-16#/25#	Φ6.4-17#	17#/24#-ID8*1
Φ4.8	15#/25#	Φ4.8-Φ6.4	15#/25#-17#/24#	Φ6.4-24#	24#-ID6*1
Φ6.4	17#/24#	Φ6.4-Φ7.9	24#-35#/36#	Φ7.9-35#	35#/36#/73#-ID8*1
Φ7.9	18#/35#	Φ7.9-Φ6.4	35#-24#	Hard tube connector	Compatible tube
Φ9.6	36#/73#	Φ9.6-Φ6.4	36#-24#	Φ1.6-14#-A	1*1/2*1/14#-OD 1/16 "
Φ12.7	82#	Φ4.8-Φ7.9	15#/25#-18#/35#	Φ3.2-16#-A	3*1/4*1/16#/25#-OD 1/8 "
T-shaped connector	Compatible tube	Φ4.8-Φ9.5	15#/25#-36#/73#	Φ1.6-Φ3.2-A	1*1/2*1/14#-OD 1/8 "
Φ3.2	16#	Φ6.4-Φ12.7	17#/24#-82#	Right-angle connector	Compatible tube
Φ4.8	15#/25#	Φ9.5-Φ12.7	36#/73#-82#	Φ6.4-24#	24#-24#
Φ6.4	17#/24#	Reducer threaded connector	Compatible tube	Φ7.9-35#	35#/36#-35#/36#
Φ7.9	18#/35#	Φ4.8-Φ3.2	15#/25#-ID4*1	-	-
Φ12.7	82#	-	-	-	-

