


Application

- ⊙ Experimental Research
- ⊙ Household
- ⊙ Detergent Conveying
- ⊙ Ink Conveying



Pump head



Features

- Compatible with various OEM peristaltic pump heads, with multiple flow options
- Easy operation, stable and reliable performance
- One-key knob operation for on/off control and speed adjustment
- Comes with a foot pad for placement on a desktop, can also be wall-mounted with optional installation plate to save space

Functions

Control Method	Potentiometer knob
Speed Resolution	0.1rpm
Display	/
Program Function	/
External Control Interface	PWM controller (50%-100% speed)

Technical Specifications

Pump Head	104KA	MN3	MN4	503Z	603K	TH504	203KA	803Z	353K
Drive Dimension	127*84.7*85mm		142*134*100mm				217*159*130mm		
Net Weight	365g	263g	830g	895g	925g	820g	980g	1490g	1840g
Flow Rate (ml/min)	≤260	≤555	≤1400	≤2185	≤2940	≤1560	≤470	≤5340	≤3915
Power Supply	DC24V,2A		DC24V,2A				DC24V,3.5A		
Power	≤24W		≤48W				≤48W		
Work Environment	Temperature range: 0°C (no freezing) to 60°C, relative humidity <85%RH								
Output Pressure	≤0.17MPa								
Motor Type	DC reducer motor with PWM speed controller								
Shell Material	ABS plastic								
Ingress Protection	IP31								

MN4 pump head



503Z pump head



353K pump head

