


Applications

- ◎ Analysis equipment
- ◎ Environmental protection
- ◎ Automatic nickel addition device
- ◎ Fermentation tank
- ◎ Cell culture



Pump Size

TH408/BT




Features

- Metal pump head, compact and exquisite in size
- Adjustable pressure block, convenient and quick tube installation
- Over-pressure gap adjustable with spring self-adaptive function
- Single pump head can accommodate 1 or 2 tubes to form single or double channels.
- Rotor set with options of 4 or 8 rollers, multiple flow rate options
- Compatible with 1.6mm wall thickness tubing

Technical Specifications

Pump Model	TH408/BT
Motor Type	Stepper
Power Supply	DC 24V
Power	≤12W
Net Weight	920g
Work Environment	Temperature range: 0℃ (no freezing) to 60℃, relative humidity <85%RH
Output Pressure	≤0.17MPa
4 Rollers – Flow Rate (unit:ml/min)	
rpm Tube#	300r
14#	45
16#	160
19#	95
25#	320

Single Channel



Double Channel

